


# Smartphone Stand Instructions

---


There are several labs in this course that ask you to record videos of the experiment to use for motion tracking. The following instructions show how to use the supplies in your lab kit to make a simple smartphone stand for filming.

---

Gather supplies.




Wrap 2 rubber bands around the wood block.



### Vertical Video Configuration

Insert the phone between the rubber bands and the block. Add the steel hex nut or other weight if the phone tips over or needs additional stabilization.



### Horizontal Video Configuration

Insert the phone between the rubber bands and the block. Add the steel hex nut or other weight if the phone tips over or needs additional stabilization.

